

Surgirack®

Efficiency, stability and transparency.

Up to 30% of laparoscopic procedural time is wasted on inefficient instrument exchanges.

Instrument exchanges are a major source of flow disruption and flow disruption significantly increases surgical errors. Surgirack uniquely arranges instruments in an ergonomic method. This greatly reduces time spent retrieving instruments. By reducing instrument exchange time, overall procedural time and surgeon frustration are reduced allowing for better surgical outcomes.

Cost Advantages: In addition to reducing procedural times, surgical instruments are protected and damage is reduced. Dislodgement of the Surgirack is prevented by 2-point conformed fixation. No more replacing expensive energy devices that fall to the floor.

In addition, the Surgirack is designed to the correct length preventing damage to fragile energy device tips which is another reason for replacement of an expensive energy device during a procedure.

Sustainability: The reusable body is autoclavable, Sterrad and EO compatible reducing instrument holder waste by 70%. A sterile, transparent section is used for each procedure.


Contact: Reach out to your local Surgirack distributor or contact JB Surgical directly via website surgirack.com.au

Dependable 2 point fixation you can rely on ▶

◀ Fast and predictable instrument retrieval. Instruments remain separated at all times

◀ Surgirack is manufactured in Australia by an Australian owned and operated company. Surgirack® is a registered trademark

◀ Available for both general surgery and bariatric length instruments



▶ Diverting instrument cords prevents entanglement

▶ Instrument jaws are always visible

▶ Curved to patient's anatomy for stability

